



TOP TRAY




BOTTOM TRAY

Flat Spanner




I-SP-X (system specific)

Depth Gauge



I-DG-24

Torque Wrench



I-TWS

Torque Attachment



I-TWS-B45  
I-TWS-B100

Wrench Converters



I-WI-CST  
For Handpiece inserts (latch-type) featuring the W&H hex.



I-WI-SL  
For Handpiece inserts (latch-type) without the W&H hex.



I-WI-SS  
For SQUARE connection of fixture mounts and instruments.



I-WI-SH  
For HEX connection of fixture mounts.